

# Koc University

## Medical Image Analysis Homework 1

Alp Akkanlar 76022  
Furkan Tuna 69730

### Part1 : ObtainForegroundMask

#### Pseudocode:

Function obtain\_foreground\_mask(image\_path)

Load image, convert to grayscale, apply Gaussian blur and Otsu's thresholding, return mask.

Function evaluate\_foreground\_mask(estimated\_mask, gold\_standard\_mask)

Calculate and return precision, recall, and F-score based on true positives, false positives, and false negatives.

Main

For each image from 1 to 3, perform the following:

Obtain foreground mask.

Load gold standard mask.

Evaluate the mask and print metrics.

Save the evaluated mask.

Handle exceptions if any occur.

#### List of Parameters:

-Threshold -> Chosen through empirical testing and analysis of the image histograms.

Otsu's method was used for its adaptive calculation of the threshold, which is particularly effective when there's a clear bi modal distribution of pixel intensities.

-Gaussian Blur: The kernel size was selected based on the level of noise present in the images. A 5x5 kernel was found to be effective in smoothing without overly blurring cell boundaries.

-Morphological Operations: The structure and size of the kernel used for operations like opening and closing, which can refine the mask by removing noise and filling holes. The kernel sizes for opening and closing were set to remove small artifacts (such as specks of noise) and close small holes within cells, respectively.

These sizes were determined by the approximate size of noise artifacts and gaps in the cell structures observed during initial inspections.

#### Table of Quantitative Metrics:

Image	Precision	Recall	F-score
Image 1	0.85	0.90	0.875
Image 2	0.80	0.88	0.840
Image 3	0.78	0.85	0.815

## Visual Results:

Image 1:




Image 2:





Image 3:



## Part 2: FindCellLocations

### Pseudocode:

Function read\_gold\_standard(gold\_path)

Load and return gold standard labels from file.

Function detect\_cell\_boundaries(rgb\_image)

Convert to grayscale, detect edges using Canny.

Function calculate\_distance\_transform(foreground\_mask, edges)

Combine mask and edges, compute and return distance transform.

Function find\_cell\_locations(dist\_transform)

Identify and return regional maxima coordinates in distance transform.

Function calculate\_metrics(detected\_labels, gold\_labels)  
Compare detected and gold labels to calculate and return precision, recall, and F-score.

Function process\_image(image\_path, mask\_path, gold\_standard\_path)  
Load image, foreground mask, and gold standard.  
Compute edges, distance transform, and detect cell locations.  
Map detected locations back to cell labels.  
Calculate metrics, print and save results.

Main  
Process each image with associated mask and gold standard, display and save output.

List of Parameters:

1. image\_path: Path to the RGB image file.
  2. mask\_path: Path to the binary foreground mask file.
  3. gold\_standard\_path: Path to the file containing gold standard cell labels.
  4. threshold in detect\_cell\_boundaries(): Edge detection threshold for Canny (default values 50 and 150).
  5. min\_distance in find\_cell\_locations(): Minimum distance between regional maxima (default 20 pixels).
  6. threshold\_abs in find\_cell\_locations(): Absolute intensity threshold for regional maxima detection (default 2).
- Edge Detection Thresholds: The Canny thresholds were chosen to effectively identify sharp transitions in grayscale intensity, which are indicative of cell boundaries. The lower threshold is relatively low to ensure weaker edges aren't missed, while the higher threshold prevents the detection of insignificant features.
  - Distance Transform Calculation: The bitwise operation combines the foreground mask and edge detection results to ensure that the distance transform calculates distances from the edges of cells, enhancing the detection of central cell points.
  - Regional Maxima Parameters: min\_distance ensures that detected cell centers are spaced out enough to correspond to distinct cells, preventing multiple detections within a single cell. The threshold\_abs helps to ignore trivial maxima in the distance transform, focusing on significant peaks that more likely represent cell centers.

Table of Quantitative Metrics:

Image	Precision	Recall	F-score
im1	1.0000	0.3471	0.5153
im2	1.0000	0.4084	0.5799
im3	1.0000	0.3615	0.5310

## Visual Results

Image1:




Image 2:





Image3

## Part3: FindCellBoundaries

### Pseudocode:

Function read\_gold\_standard(gold\_path)

    Load gold standard data from a file located at gold\_path.

Function detect\_cell\_boundaries(rgb\_image)

    Convert RGB image to grayscale and apply the Canny algorithm to detect edges.

Function calculate\_distance\_transform(foreground\_mask, edges)

    Combine foreground mask and edge detection results using a bitwise OR, then compute and return the distance transform.

Function find\_cell\_locations(dist\_transform)

Identify and return the coordinates of regional maxima in the distance transform as potential cell locations. -> Because we didn't use jupyter notebook we had to save the cell locations in a json format in order to get them to another script.

Function calculate\_metrics(detected\_labels, gold\_labels)

Calculate precision, recall, and F-score by comparing detected labels with gold standard labels.

Function process\_image(image\_path, mask\_path, gold\_standard\_path)

Load the image and mask, process to find edges and calculate the distance transform.

Detect cell locations, calculate and evaluate metrics against the gold standard.

Output visual results and metrics.

Main

For each image, mask, and gold standard pair, invoke process\_image to perform detection and evaluation.

## List of Parameters:

1. gold\_path:
    - Purpose: Path to the file containing gold standard labels.
    - Selection: This should be the complete path to ensure the program can correctly locate and load the necessary file.
  2. rgb\_image:
    - Purpose: The RGB image on which cell boundary detection will be performed.
    - Selection: Passed directly to the function to apply image processing techniques specific to the image's characteristics.
  3. foreground\_mask:
    - Purpose: A binary mask indicating foreground (cell) areas.
    - Selection: Typically derived from a previous step that isolates the cells from the background, affecting subsequent distance transform calculations.
  4. edges:
    - Purpose: The edge-detected output from the Canny algorithm.
    - Selection: Canny edges help refine the distance transform by incorporating structural information from the cells.
  5. dist\_transform:
    - Purpose: Distance transform of the combined foreground mask and edges.
    - Selection: This transform helps identify regional maxima as potential cell centers based on their distance from the nearest edge, critical for accurate cell location identification.
  6. detected\_labels and gold\_labels:
    - Purpose: Arrays of detected and gold standard labels for metric calculation.
    - Selection: Used to assess the accuracy of the segmentation against a verified standard, influencing how the detection algorithm's performance is quantified.
- Edge Detection Thresholds (low\_threshold, high\_threshold for Canny):
    - These thresholds are crucial for effective edge detection. Low threshold values may result in detecting too much noise, whereas too high values might miss significant edges. Typically, a ratio of 1:2 or 1:3 between the lower and upper thresholds is effective.
  - Distance Transform Method (cv2.DIST\_L2):
    - The choice of using the L2 norm (Euclidean distance) for distance transforms provides a more natural measurement of distance as opposed to other methods like L1 (Manhattan distance), especially relevant in circular or spherical cell structures where Euclidean distance gives a realistic portrayal of cell boundaries.
  - Regional Maxima Detection Parameters (min\_distance, threshold\_abs):
    - min\_distance: This controls the minimum distance between detected centers, preventing too closely situated maxima that might correspond to the same cell. It should be set based on the expected minimum cell size.

- `threshold_abs`: This threshold helps ignore smaller maxima that might be due to noise or non-cell elements, ensuring only significant maxima are considered.
- Seed Points for Region Growing:
  - Seed points are chosen based on regional maxima of the distance transform, which are expected to be central within cells. Ensuring these are well-selected directly impacts the accuracy and efficiency of the region-growing algorithm.

## Quantitative Metrics:

Image	Precision	Recall	F1 Score	IOU (0.5)	Dice (0.5)	IOU (0.75)	Dice (0.75)	IOU (0.9)	Dice (0.9)
im1	0.79	0.05	0.09	0.64	0.00	0.00	0.00	0.00	0.00
im2	0.77	0.08	0.14	0.00	0.00	0.00	0.00	0.00	0.00
im3	0.65	0.07	0.13	0.00	0.00	0.00	0.00	0.00	0.00

Our segmentation algorithm's performance metrics indicate that while the cells it identifies are generally correct (reflected by relatively high precision), it's only recognizing a small fraction of the actual cells present (hence the low recall). This suggests that the algorithm is overly conservative in its identification process. The F1 Score, which balances precision and recall, is low, signifying that the algorithm's cautious approach is causing it to miss too many true cell identifications, thus reducing its overall effectiveness. When we consider the IoU and Dice metrics, particularly at higher thresholds, the scores drop to zero, indicating that the algorithm's accuracy diminishes significantly as the criteria for a match become stricter. This could be a result of the algorithm not capturing the full extent of each cell, possibly due to variations in cell shape and size that are not accommodated by the current segmentation parameters.

## Visual Results:

Image 1:




Image 2:





Image 3:



## Part 4: Segmentation of blood vessels in fundus photography images

### Pseudocode:

```
Function preprocess_image(image_path)
    Load image in grayscale
    Apply CLAHE (Contrast Limited Adaptive Histogram Equalization) to the image
    Return preprocessed image
```

```
Function segment_vessels(image)
```

```

Apply Gaussian and median blurring to smooth the image
Compute Sobel gradients in x and y directions
Calculate gradient magnitude
Apply local thresholding to get vessel mask
Clean the mask by removing small objects and holes
Return the cleaned vessel mask

```

```

Function calculate_metrics(estimated_mask, gold_mask)
    Convert masks to binary format
    Calculate true positives (TP), false positives (FP), and false negatives (FN)
    Compute precision, recall, and F1 score from TP, FP, FN
    Return precision, recall, F1 score

```

```

Function process_images(image_paths, gold_paths)
    For each image and corresponding gold standard path
        Preprocess the image
        Segment vessels from the preprocessed image
        Load gold standard mask
        Print unique values and dimensions of gold and estimated masks
        Calculate and print precision, recall, and F1 score
        Save segmented image

```

```

Main
    Define image paths and corresponding gold standard paths
    Call process_images with image and gold paths

```

### List of Parameters:

1. CLAHE parameters: clipLimit and tileGridSize determine the contrast enhancement limit and the number of tiles in the grid for histogram equalization. The values are chosen based on empirical results, optimizing the visibility of vessels while minimizing noise.
2. Gaussian and Median Blur: The kernel sizes for blurring filters are chosen to sufficiently smooth the image while preserving the vessel structures.
3. Sobel Operator: The kernel size for the Sobel operator impacts the edge detection sensitivity. A kernel size of 5 is a balance between finding fine details and avoiding noise.
4. Local Thresholding: block\_size in local thresholding affects the adaptiveness of the threshold. The size is selected based on the typical size of the vessels to ensure that larger areas are not mistakenly thresholded as vessels.
5. Cleaning Parameters: min\_size for removing small objects and holes is based on the approximate size of the noise artifacts to be removed.


### Quantitative Metrics:

Image	Precision	Recall	F1 Score
d4_h.jpg	0.2159	0.5980	0.3173
d7_dr.jpg	0.1596	0.4959	0.2414
d11_g.jpg	0.1516	0.5284	0.2356




**Visual Results:**

**D4\_h segmented**



**D7\_dr segmented**



**D11\_g segmented**

